



Hewlett Packard
Enterprise

HPE Security ArcSight Compliance Insight Package for HIPAA 3.0

Software Version: ESM 6.8c, 6.9.1c, or 6.11.0

Release Notes

November 20, 2017

Legal Notices

Warranty

The only warranties for Hewlett Packard Enterprise products and services are set forth in the express warranty statements accompanying such products and services. Nothing herein should be construed as constituting an additional warranty. Hewlett Packard Enterprise shall not be liable for technical or editorial errors or omissions contained herein.

The information contained herein is subject to change without notice.

The network information used in the examples in this document (including IP addresses and hostnames) is for illustration purposes only.

HPE Security ArcSight products are highly flexible and function as you configure them. The accessibility, integrity, and confidentiality of your data is your responsibility. Implement a comprehensive security strategy and follow good security practices.

This document is confidential.

Restricted Rights Legend

Confidential computer software. Valid license from Hewlett Packard Enterprise required for possession, use or copying. Consistent with FAR 12.211 and 12.212, Commercial Computer Software, Computer Software Documentation, and Technical Data for Commercial Items are licensed to the U.S. Government under vendor's standard commercial license.

Copyright Notice

© Copyright 2017 Hewlett Packard Enterprise Development, LP

Follow this link to see a complete statement of copyrights and acknowledgements:

<https://community.saas.hpe.com/t5/Discussions/Third-Party-Copyright-Notices-and-License-Terms/td-p/1589228>

Support

Contact Information

Phone	A list of phone numbers is available on the HPE Security ArcSight Technical Support Page: https://softwaresupport.hpe.com/support-contact-information
Support Web Site	https://softwaresupport.hpe.com
Protect 724 Community	https://community.saas.hpe.com/t5/ArcSight/ct-p/arc sight

Contents

- CIP for HIPAA 3.0 4
 - Release Contents 4
 - Installing CIP for HIPAA 3.0 4
 - Performance Impact of CIP for HIPAA 5
 - Open Issue in this Release 5

- Send Documentation Feedback 6

CIP for HIPAA 3.0

In 1996, the Health Insurance Portability and Accountability Act (HIPAA) was enacted by the United States Congress. Title II of HIPAA has standardized the way electronic health care information is transmitted among providers and insurers by defining a set of standard transactions. In addition, it defines a Security Rule which specifies that covered entities (such as health insurers, hospitals, medical providers) put appropriate safeguards in place to protect the Electronic Protected Healthcare Information (EPHI). The Compliance Insight Package for HIPAA coupled with ArcSight ESM can assist you in complying with the HIPAA CIP Standards.

Release Contents

The files included in this release are:

File name	Description
ESM_HIPAA_Solution_ReleaseNotes_3.0.pdf	This document. Product descriptions and open issues
ESM_HIPAA_SolutionGuide_3.0.pdf	<i>ArcSight ESM Compliance Insight Package for HIPAA 3.0 Solution Guide</i> Product architecture, installation, configuration, and operation instructions with a description of product contents.
ComplianceInsightPackage-HIPAA.3.00049.arb	Installable package bundle for all operating systems. Contains all the resources of the Compliance Insight Package for HIPAA. Note: If you use Internet Explorer to download the ARB file, it may convert the ARB file to a ZIP file. If this occurs, rename the ZIP file back to an ARB file before importing into ArcSight ESM.

Installing CIP for HIPAA 3.0

The Compliance Insight Package for HIPAA (CIP for HIPAA) is supported on the following versions of ESM:

- ESM 6.8c
- ESM 6.9.1c
- ESM 6.11.0

CIP for HIPAA 3.0 is for fresh installation only.

Complete installation and configuration instructions for the Compliance Insight Package for HIPAA are located in the document called *ArcSight ESM Solution Guide Compliance Insight Package HIPAA 3.0* (ESM_HIPAA_SolutionGuide_3.0.pdf).

Performance Impact of CIP for HIPAA

ArcSight ESM Solution Packages contain data monitors and rules that can place an additional load on the ArcSight ESM Manager, which may impact the ArcSight ESM Manager performance. If your ESM system is operating at an average event per second (EPS) rate that has maximized the CPU utilization, you might experience a reduced average EPS rate after installing the CIP for HIPAA package. If this performance impact occurs, ArcSight recommends that you disable unneeded data monitors and rules to reduce the load on the Manager.

Open Issue in this Release

The following issue is open in this release.

Number	Description and Workaround
SOL-5106	On high EPS rates, the "Successful Attack - Brute Force Login" rule could affect ESM Manager performance. Workaround: If this performance impact occurs, ArcSight recommends that you disable unneeded data monitors and rules to reduce the load on the ESM Manager.

Send Documentation Feedback

If you have comments about this document, you can [contact the documentation team](#) by email. If an email client is configured on this system, click the link above and an email window opens with the following information in the subject line:

Feedback on Release Notes (Compliance Insight Package for HIPAA 3.0 ESM 6.8c, 6.9.1c, or 6.11.0)

Just add your feedback to the email and click send.

If no email client is available, copy the information above to a new message in a web mail client, and send your feedback to arc-doc@hpe.com.

We appreciate your feedback!